

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

P85a; PIK3R1 Fusion Protein



Basic Information

Catalog Number:

Ag2344

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RDTADGTFLVRDASTKMHGDYTLTRKGGNNKLIKIFH
RDGKYGFSDPLTFSSVVELINHRYNESLAQYNPKLDVK
LLYPVSKYQQDQVVKEDNIEAVGKKLHKYNTQFQEK
SREYDRLYEYTRTSQEIQMKRTAIEAFNETIKIFEEQC
QTQERYSKYIEKFKREGNEKEIQRMHNYDKLSRISE
IIDSRRLEEDLKKQAAEYREIDKRMNSIKPDLIQLRKR
DQYLMWLTQKGVQRQKLNELGNENTEDQYSLVED
DEDLPHHDEKTWNVGSSNRNKAENLLRGKRDGTFLVR
ESSKQGCYACSVVVDGEVVKHCVINKTATGYGFAEPY
NLYSSLKELVLHYQHTSLVQHNDLSLNTLAYPVYAQQ
RR

(78-454 aa encoded by BC030815)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.