

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CUL2 Fusion Protein



Basic Information

Catalog Number:

Ag28654

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

ESVLEVHGKVFVQLINTVLNGDQHFMSALDKALTSVV
NYREPKSVCKAPELLAKYCDNLLKKSAGMTENEVED
RLTSFITVFKYIDDKDVFQKFYARMLAKRLIHGLSMSM
DSEEMINLKLQACGYEFTSKLHRMYTDM SVSADLN
NKFNNFIKNQDTVIDLGISFQIYVLQAGAWPLTQAPS
STFAIPQELEKSVQMFELFYSQHFSGRKLTWLHLYLCT
GEVKMNYLGKPYVAMVTTYQMAVLLAFNNSETVSYK
ELQDSTQMNEKELTKTIKSLLDVKMINHDSEKEDIDAE
SSFSLNMNFSSKRTKFKITTSMQKDTPEMEQTRSAV
DEDKMYLQAAIVRIMKARKVLRHNALIQEVISQSRA
RFNPSJSMIKKCI EVLIDKQYIERSQASADEYSYVA
(341-745 aa encoded by BC009591)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.