

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Snf1lk Fusion Protein



## Basic Information

Catalog Number:

Ag3913

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LDTNCSGVFRHRSISPSLLDTAISEEARQGPSLEEEQ  
EVQEPPLPGSTGRRHTLAEVSTHFSPLNPPCIIVSSSAT  
ASPSEGTSSDSCLPFSASEGPAGLGSGLATPGLLTSS  
PVRLASPFSGQSATPVLTQAGLGTAVLPPVSFQEG  
RRASDTSLTQGLKAFRQQLRKNARTKGFGLNKKGLA  
RQVCQSSVRTPRGGMSTFHTPAPSSGLQGCTTSNRE  
GRSLLEEV LHQQRLLQLQHHSSTAAA SSGCQ QGPQL  
SPVPYVLAPCDLLVSGIPLLPTLLQAGMSPVASAAH  
LLDTHLHISAGPVALPTGPLPQCLTRLSPGCDPAGLPQ  
GDCEMEDLTS

(423-769 aa encoded by NM\_010831)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.