

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

ACK1 Fusion Protein



Basic Information

Catalog Number:

Ag5419

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MQPEEGTGWLLELLSEVQLQQYFLRLRDDLNVTRLSH
FEYVKNEDLEKIGMRPGQRRLWEAVKRRKALCKRKS
WMSKVFSGKRLEAEFPPHHSQSTFRKTSPAPGGPAG
EGPLQSLTCLIGEKDLRLEKLGDSFGVVRGEWDAP
SGKTVSVAVKCLKPDVLSQPEAMDDFIREVNAMHSL
DHRNLIIRLYGVVLTTPMKMVTTELAPLGSLLDRLRKHQG
HFLLGTLTRYAVQVAEGMGYLESKRFIHRDLAARNLLL
ATRDLVKIGDFGLMRALPQNDDHYVMQEHKRVFPA
WCAPESLKRTRTFSHASDTWMFGVTLWEMFTYGGQP
WIGLNGSQILHKIDKEGERLPRPEDCPQDIYNVMVQC
WAHKPEDRPTFVALRDFL

(1-385 aa encoded by BC028164)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.