

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# cTnT Fusion Protein



## Basic Information

Catalog Number:

Ag7716

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSDIEEVVEEYEEEQEEAAVEEQEEAAEEDAEAEA  
ETEETRAEEDDEEEEAKEAEDGPMEESEKPKPRSFMPNL  
VPPKIPDGERVDFDDIHRKRMEKDLNELQALIEAHFEN  
RKKEEEELVSLKDRIERRRAERAEOQRIRNEREKERQN  
RLAEERARRDEEENRRKAEDEARKKKALSMMHFGGY  
IQKTERKSGKRQTEREKKKKILAERRKVLAIIDLHNLNEDQL  
REKAKELWQSIYNLEAEKFDLQEKFKQQKYEINVLNRN  
INDNQKVSCTRGGAKVTGRWK  
(1-285 aa encoded by BC002653)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.