

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

ABI3 Fusion Protein



Basic Information

Catalog Number:

Ag8164

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAELQQLQEF E I P T G R E A L R G N H S A L L R V A D Y C E D N Y
V Q A T D K R K A L E E T M A F T T Q A L A S V A Y Q V G N L A G H T L
R M L D L Q G A A L R Q V E A R V S T L G Q M V N M H M E K V A R R E
I G T L A T V Q R L P P G Q K V I A P E N L P L T P Y C R R P L N F G C L
D D I G H G I K D L S T Q L S R T G T L S R K S I K A P A T P A S A T L G R P
P R I P E P V H L P V V P D G R L S A A S S A S L A S A G S A E G V G G
A P T P K G Q A A P P A P L P S S L D P P P P P A A V E V F Q R P P T L
E E L S P P P D E E L P L D L P P P P L D G D E L G L P P P P G F G P
D E P S W W P A S Y L E K V V T L Y P Y T S Q K D N E L S F S E G T V I C
V T R R Y S D G W C E G V S S E G T G F F P G N Y V E P S C
(1-366 aa encoded by BC007780)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.