

For Research Use Only

Goat anti-Rabbit IgG

Catalog Number: B900210



www.ptgcn.com

Physical State	Liquid
Protein Concentration	1 mg/ml
Size	100ug
Purification	The antibody was purified from antisera by immunofinity chromatography.
Buffer	PBS, 50% glycerol, pH 7.3
Note	This product is for research use only, not for diagnostic or therapeutic use.
Preservative	0.02% Sodium Azide
Storage	-20°C, avoid repeated thawing and freezing

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.