

For Research Use Only

CoraLite® Plus 488-conjugated RPLPO Monoclonal antibody



Catalog Number: CL488-68056

Basic Information

Catalog Number:

CL488-68056

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1829

GenBank Accession Number:

BC009867

GeneID (NCBI):

6175

UNIPROT ID:

P05388

Full Name:

ribosomal protein, large, PO

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Protein G purification

CloneNo.:

2A2E1

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human, Mouse, Rat

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human RPLP0 (CL488-68056, Clone:2A2E1) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).