

For Research Use Only

CoraLite®594-conjugated DMC1 Monoclonal antibody



Catalog Number:CL594-67176

Basic Information

Catalog Number:

CL594-67176

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6951

GenBank Accession Number:

BC035658

GeneID (NCBI):

11144

UNIPROT ID:

Q14565

Full Name:

DMC1 dosage suppressor of mck1
homolog, meiosis-specific
homologous recombination (yeast)

Calculated MW:

165 aa, 18 kDa, 38 kDa

Observed MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

2D5C9

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

Human, mouse, rat

Positive Controls:

IF : mouse testis tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite®594 DMC1 antibody (CL594-67176, Clone: 2D5C9) at dilution of 1:200.