

For Research Use Only

# AFP Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP00001-2

## Capture Antibody Information

Catalog Number:  
82982-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230435G1

Reactivity:  
Human

GenBank:  
BC027881

Conjugate:  
Unconjugated

Full name:  
alpha-fetoprotein

Gene ID:  
174

## Detection Antibody Information

Catalog Number:  
82982-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230435F11

Reactivity:  
Human

GenBank:  
BC027881

Conjugate:  
Unconjugated

Full name:  
alpha-fetoprotein

Gene ID:  
174

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00001-2 targets AFP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AFP Recombinant antibody, PBS Only (Capture) 82982-3-PBS (230435G1). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: AFP Recombinant antibody, PBS Only (Detector) 82982-2-PBS (230435F11). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative AFP matched antibody pairs: MP00001-1, MP00001-3, MP50072-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00001-2, AFP Recombinant Matched Antibody Pair - BSA and Azide free, measured by cytometric bead array. Capture antibody: 82982-3-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 1.25-20 ng/mL