

For Research Use Only

NLRP1 Polyclonal antibody

Catalog Number: 24792-1-AP



Basic Information

Catalog Number: 24792-1-AP	GenBank Accession Number: BC051787	Purification Method: Antigen affinity purification
Size: 800 µg/ml	GeneID (NCBI): 22861	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q9C000	
Isotype: IgG	Full Name: NLR family, pyrin domain containing 1	
Immunogen Catalog Number: AG20516	Calculated MW: 1473 aa, 166 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human tonsillitis tissue, human lymphoma tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 24792-1-AP (NALP1 antibody at dilution of 1:50) (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 24792-1-AP (NALP1 antibody at dilution of 1:50) (under 40x lens).