

For Research Use Only

# V5-tag Monoclonal antibody

Catalog Number: 66007-1-Ig **37 Publications**



## Basic Information

<b>Catalog Number:</b> 66007-1-Ig	<b>GenBank Accession Number:</b> GeneID (NCBI): 4	<b>Purification Method:</b> Protein G purification
<b>Concentration:</b> 800 µg/ml	<b>Full Name:</b> V5 Tag	<b>CloneNo.:</b> 2B7A9
<b>Source:</b> Mouse	<b>Calculated MW:</b> 1 kDa	
<b>Isotype:</b> IgG1		

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB, IP, CoIP, RIP

**Species Specificity:**  
recombinant protein

## Background Information

The antibody is developed against a synthetic 14 aa peptide (GKPIPPLLGLDST), the V5 tag, which presents in many vectors.

## Notable Publications

Author	Pubmed ID	Journal	Application
Jinbo Fan	32967926	Life Sci Alliance	WB
Xing Sun	30237395	Cell Death Dis	
Ke Zheng	31777585	J Cancer	

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data