

For Research Use Only

ACF1 Monoclonal antibody

Catalog Number: 67947-1-Ig



Basic Information

Catalog Number:

67947-1-Ig

Size:

1000 μ g/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30607

GenBank Accession Number:

NM_013448

GeneID (NCBI):

11177

UNIPROT ID:

Q9NRL2

Full Name:

bromodomain adjacent to zinc finger domain, 1A

Calculated MW:

179 kDa

Purification Method:

Protein A purification

CloneNo.:

1C4F8

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, HEK-293 cells, Jurkat cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67947-1-Ig (ACF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.