

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

YAP1 Fusion Protein

Basic Information

Catalog Number:

Ag28194

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6[×]His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Note: Remove endotoxin and sterilization is required for cell assay.

Validated Application:

Blocking peptide

Peptide Sequence:

PTAQLRQSSFEPDDVPLPAGWEMAKTSSGQRYFL
NHIDQTTTWQDPRKAMLSQMNVTAPTSPVQQNM
MNSASGPLPDGWEQAMTQDGEIYYINHKNKTTSWLD
PRLDPRFAMNQRISSQAPVKQPPPLAPQSPQGGVMG
GSNSNQQQMRLQQLQMEKERLRLKQQLLRQAMR
NINPSTANSPKQCEALRSLQPTLEQDGGTQNPVSSP
GMSQELRMTTNSDDPFLNSGTYHSRDESTDSGLSMS
SYSVPRTPDDFLNSVDEMDTGDTINQSTLPSQQNRF
DYLEAIPGTNVDLGTLEGDMNIEGELMPSLQEQALS
DILNDMESVLAATKLDKESFTWL
(155-504 aa encoded by BC038235)

Reconstitution and Storage

Reconstitution:

Reconstitute at 0.25 µg/ µl in 200 µl sterile water for short-term storage.

Reconstitution with 200 µl 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details).

If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used).

Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

Shipping:

The product is shipped at ambient temperature.

Upon receipt, store it immediately at the recommended temperature (see below).

Stability and Storage

Store for up to 12 months at -20°C to -80°C as lyophilized powder.

Storage of Reconstituted Protein


Short Term Storage:

Store at 2-8°C for (1-2 weeks).

Long Term Storage:

Aliquot and store at -20°C to -80°C for up to 3 months, buffer containing 50% glycerol is recommended for reconstitution. Avoid repeat freeze-thaw cycles.

Selected Validation Data



For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.