

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human A1BG protein



## Basic Information

Catalog Number:

Ag5414

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MAPVSWITPGLKTTAVCRGVLRGVTFLLRREGDHEFL  
EVPEAQEDVEATFPVHQPGNYSCSYRTDGEGALSEP  
SATVTIEELAAPPPVLMHHGESSQVLHPGNKVTLC  
VAPLSGVDFQLRRGKELLVPRSSTSPDRIFFFHLNAVA  
LGDGGHYTCRYRLHDNQNGWSGDSAPVELLSDETLP  
APEFSPPEPESGRALRLRCLAPLEGARFALVREDRGGRR  
VHRFQSPAGTEALFELHNISVADSANYSCVYVDLKPP  
FGGSAPSERLELHVDGPPRPQLRATWSGAVLAGRD  
AVLRCEGPIPDVTFELLREGETKAVKTVRTPGAAANLE  
LIFVGPQHAGNYRCRYRSWVPHTFESESDPVELLVA  
ES

(1-373 aa encoded by BC035719)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.