

For Research Use Only

# CoraLite® Plus 488 Anti-Human CD3 (OKT3)



Catalog Number: **CL488-65133**

## Basic Information

<b>Catalog Number:</b> CL488-65133	<b>GenBank Accession Number:</b> BC049847	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100tests, 5 µl/test	<b>GeneID (NCBI):</b> 916	<b>CloneNo.:</b> OKT3
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000198851	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Isotype:</b> IgG2a, kappa	<b>UNIPROT ID:</b> P07766	
	<b>Full Name:</b> CD3e molecule, epsilon (CD3-TCR complex)	
	<b>Calculated MW:</b> 207 aa, 23 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

CD3 is a multimeric protein associated with the T-cell receptor (TCR) to form a complex involved in antigen recognition and signal transduction (PMID: 15885124). CD3 is composed of CD3  $\gamma$ ,  $\delta$ ,  $\epsilon$ , and  $\zeta$  chains (PMID: 1826255). It is expressed by thymocytes in a developmentally regulated manner, T cells, and some NK cells (PMID: 3289580). The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction (PMID: 11985657).

## Storage

**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with APC-Anti-Human CD19 (APC-65145, Clone: SJ25C1) and 5.00 ul CoraLite® Plus 488-conjugated Anti-Human CD3 (CL488-65133, Clone: OKT3). Cells were not fixed.



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with APC-Anti-Human CD19 (APC-65145, Clone: SJ25C1) and CoraLite® Plus 488-conjugated Mouse IgG2a isotype control antibody. Cells were not fixed.