

For Research Use Only.
Not For Use In Diagnostics.

CoraLite594 – conjugated Goat Anti-Mouse IgG(H+L)

Catalog Number:SA00013-3



www.ptglab.com

Information

Catalog Number:

SA00013-3

Host:

Goat

Applications:

IF, FC

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

CoraLite®594

Recommended Dilutions

1:100 – 1:500

Preservative

0.02% Sodium Azide

Fluorophore

CL594 Fluorescent Dye
Amax=593 nm, Emax=614 nm

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using Lamin B1 antibody (66095-1-Ig, Clone: 3C10G12) at dilution of 1:400 and SA00013-3 CoraLite®594-Conjugated Goat Anti-Mouse IgG(H+L).