


p53 SIGNALING

www.ptglab.com



Related Products

Legend

- Association
- Directly Activates
- Indirectly Activates
- Inhibits
- Ac Acetylation
- Ub Ubiquitin
- P Phosphorylates

Antibody Name	Catalog Number	Type	Applications
14-3-3	14503-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
ATR	19787-1-AP	Rabbit Poly	ELISA, WB
BAX	50599-2-Ig	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
BCL2	12789-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
BRC1A1	20649-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
BRC1A2	19791-1-AP	Rabbit Poly	ELISA, IHC, WB
CASP3	19677-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
CCNE1 Cyclin E1	11955-1-AP	Rabbit Poly	ELISA, FC, WB
CCNE2 Cyclin E2	11935-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
CDK4	10122-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
CDK2	11026-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
CHK1	10362-1-AP	Rabbit Poly	ELISA, FC, IF, WB
CHK2	13954-1-AP	Rabbit Poly	ELISA, IF, IP, WB

Antibody Name	Catalog Number	Type	Applications
CDKN2A	10883-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
Erk1	11257-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
Erk2	16443-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
Erk2	51068-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
GADD45A	13747-1-AP	Rabbit Poly	ELISA, IP, WB
GADD45G	10277-1-AP	Rabbit Poly	ELISA, WB
HDAC1	10197-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
HDAC2	13922-3-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
HRAS	18295-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
JNK1	51151-1-AP	Rabbit Poly	ELISA, IF, IP, WB
KAT2B/PCAF	13983-1-AP	Rabbit Poly	ELISA, IHC, IP
KRAS	12063-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
NFKB1A	51066-1-AP	Rabbit Poly	ELISA, IHC, IP, WB

Antibody Name	Catalog Number	Type	Applications
MAP2K2/MMK2	11049-1-AP	Rabbit Poly	ELISA, IF, IP, WB
p16-INKAA	10883-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
p21	10355-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
p53	10442-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
PCNA	10205-2-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
PCNA	60097-1-Ig	Mouse Mono	ELISA, FC, IF, IHC, IP, WB
PIAS3	13486-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
PRKCD	14188-1-AP	Rabbit Poly	ELISA, IP, WB
PRKDC	19983-1-AP	Rabbit Poly	ELISA, WB
Rb1	10048-2-Ig	Mouse Mono	ELISA, IHC, IP, WB
Rb1	17218-1-AP	Rabbit Poly	ELISA, WB
SIRT1	13161-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
USP5	10473-1-AP	Rabbit Poly	ELISA, IF, IP, WB

This number shows the amount of times our antibody has been cited in a publication.