


INSULIN PATHWAY

www.ptglab.com



Related Products

Legend

- Association
- Directly Activates
- > Indirectly Activates
- | Inhibits
- | Indirectly Inhibits
- Ⓟ Phosphorylates

| Antibody Name | Catalog Number | Type | Applications |
|-------------------|----------------|-------------|----------------------------|
| 4EBP1 | 60246-1-Ig | Mouse Mono | ELISA, IF, IHC, WB |
| AKT | 60203-2-Ig | Mouse Mono | ELISA, FC, IF, IHC, IP, WB |
| ATP citrate lyase | 15421-1-AP | Rabbit Poly | ELISA, FC, IF, IP, WB |
| BAD | 10435-1-AP | Rabbit Poly | ELISA, IHC, WB |
| CBLC | 12781-1-AP | Rabbit Poly | ELISA, FC, IHC, IP, WB |
| CBL | 13104-1-AP | Rabbit Poly | ELISA, WB |
| CRK | 16685-1-AP | Rabbit Poly | ELISA, IF, WB |
| EIF4B | 17917-1-AP | Rabbit Poly | ELISA, IF, IP, WB |
| EIF4E | 12227-1-AP | Rabbit Poly | ELISA, IHC, IP, WB |
| Flotillin 1 | 15571-1-AP | Rabbit Poly | ELISA, FC, IF, IP, WB |
| Flotillin 2 | 11446-1-AP | Rabbit Poly | ELISA, WB |
| GAB2 | 22549-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| GLUT4 | 21048-1-AP | Rabbit Poly | ELISA, IHC, WB |

| Antibody Name | Catalog Number | Type | Applications |
|-------------------|----------------|-------------|----------------------------|
| GRB2 | 10254-2-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| GSK3B | 22104-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| Insulin | 15848-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| IRS1 | 17509-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| JNK | 51151-1-AP | Rabbit Poly | ELISA, IF, IP, WB |
| JNK2 | 17337-1-AP | Rabbit Poly | ELISA, IF, WB |
| MTOR | 20657-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| p70 S6K | 14485-1-AP | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| PDK2 | 15647-1-AP | Rabbit Poly | ELISA, IHC, IP, WB |
| PI3K p110 (alpha) | 20583-1-AP | Rabbit Poly | ELISA, IHC, WB |
| PI3K p110 (beta) | 20584-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| PI3K p55 (gamma) | 55056-1-AP | Rabbit Poly | ELISA, IP, WB |
| PI3K p85 (alpha) | 60225-1-Ig | Mouse Mono | ELISA, IF, IHC, IP, WB |

| Antibody Name | Catalog Number | Type | Applications |
|-----------------|----------------|-------------|------------------------|
| PKC alpha | 66421-1-Ig | Mouse Mono | ELISA, IF, IHC, WB |
| PKC iota | 13883-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| PKC zeta | 26899-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| PPP1CB-Specific | 10140-2-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| PTEN | 22034-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| Raptor | 20984-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| SHP2 | 20145-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| SOCS3 | 14025-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| SOS1 | 55041-1-AP | Rabbit Poly | ELISA, IF, IP, WB |
| TC10 | 17805-1-AP | Rabbit Poly | ELISA, IF, WB |

Ⓢ This number shows the amount of times our antibody has been cited in a publication.