


HYPOXIA SIGNALING PATHWAY

www.ptglab.com



Related Products

Legend

- Activates
- ⊣ Inhibits

AntigenName	CatalogNo	Host	Application	Reactivity
4EBP1	KDNKO 13 60246-1-Ig	Mouse Mono	WB, IF, ELISA	Human
AKT	KDNKO 316 10176-2-AP	Rabbit Poly	WB, IP, IHC, IF, FC, ELISA	Human, Mouse, Rat, A. sphaerocephala + more
AKT	KDNKO 51 60203-2-Ig	Mouse Mono	WB, IP, IHC, IF, ELISA	Human, Mouse, Rat, Hamster
CL488-conjugated AKT	CL488-60203	Mouse Mono	IF, ELISA	Human, Mouse, Rat
AKT1	KDNKO 8 51077-1-AP	Rabbit Poly	WB, IHC, IF, ELISA	Human
Phospho-AKT1 (Ser473)	28731-1-AP	Rabbit Poly	WB, ELISA	Human, Mouse, Rat
AKT2	KDNKO 2 17609-1-AP	Rabbit Poly	WB, IP, IHC, ELISA	Human, Mouse, Rat
AKT3	KDNKO 7 21641-1-AP	Rabbit Poly	WB, IHC, ELISA	Human, Mouse, Rat
Phospho-AKT (Ser473)	KDNKO 179 66444-1-Ig	Mouse Mono	WB, IHC, IF, ELISA	Human, Mouse, Rat, Chicken
ARNT, HIF1B	KDNKO 7 14105-1-AP	Rabbit Poly	WB, IP, IHC, IF, ELISA	Human, Mouse, Rat
ARNT, HIF1B	66732-1-Ig	Mouse Mono	WB, IF, ELISA	Human
BCL2	KDNKO 63 60178-1-Ig	Mouse Mono	WB, IP, IHC, IF, FC, CoIP, ELISA	Human, Pig, Rat, Swine
BCL2	KDNKO 869 12789-1-AP	Rabbit Poly	WB, IP, IHC, IF, FC, chIP, ELISA	Human, Astragalus Membranaceus, + more
CBP, EIF4E	KDNKO 7 11149-1-AP	Rabbit Poly	WB, IHC, IF, FC, ELISA	Human
EIF4E	KDNKO 66655-1-Ig	Mouse Mono	WB, IHC, IF, ELISA	Human, Mouse, Rat
EIF4E2	12227-1-AP	Rabbit Poly	WB, IP, IHC, ELISA	Human, Mouse, Rat
EPO	66975-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Mouse
EPO	KDNKO 17908-14-AP	Rabbit Poly	WB, IHC, ELISA	Human
ERK1/2	KDNKO 5 51068-1-AP	Rabbit Poly	WB, IHC, IF, ELISA	Human, Mouse, Rat
ERK1/2	KDNKO 163 16443-1-AP	Rabbit Poly	WB, IHC, IF, ELISA	Human, Mouse, Rat, A. Sphaerocephala + more
Phospho-ERK1/2 (Thr202/Tyr204)	28733-1-AP	Rabbit Poly	WB, ELISA	Human, Mouse, Rat
CL594-conjugated ERK1/2	CL594-66192	Mouse Mono	IF, ELISA	Human, Mouse
CL488-conjugated ERK1/2	CL488-66192	Mouse Mono	IF, ELISA	Human, Mouse
ERK1/2	KDNKO 20 11257-1-AP	Rabbit Poly	WB, IHC, IF, ELISA	Human, Mouse, Rat
ERK1/2	KDNKO 5 67170-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Mouse, Rat
HIF1a	KDNKO 162 20960-1-AP	Rabbit Poly	WB, IP, IHC, IF, FC, CoIP, chIP, ELISA	Human, Chicken
HIF1a	KDNKO 3 66730-1-Ig	Mouse Mono	WB, IHC, IF, ELISA	Human
IKKBK	KDNKO 29 15649-1-AP	Rabbit Poly	WB, IP, IHC, IF, ELISA	Human, Mouse, Rat
IL-6	KDNKO 174 21865-1-AP	Rabbit Poly	IHC, IF, ELISA	Human, Dog, Goat
IL-6	KDNKO 15 66146-1-Ig	Mouse Mono	WB, IHC, IF, FC, ELISA	Human, Mouse
IL-6R	KDNKO 5 66855-1-Ig	Mouse Mono	WB, IHC, ELISA	Human
IL-6R	23457-1-AP	Rabbit Poly	WB, IP, IHC, ELISA	Human, Mouse, Rat

AntigenName	CatalogNo	Host	Application	Reactivity
MEK1/2	KDNKO 6 11049-1-AP	Rabbit Poly	WB, IP, IF, ELISA	Human, Mouse, Rat
MEK1-Specific	KDNKO 2 51080-1-AP	Rabbit Poly	WB, IHC, ELISA	Human
MNK1	KDNKO 2 10136-1-AP	Rabbit Poly	WB, IF, ELISA	Human, Mouse, Rat
mTOR	KDNKO 7 66888-1-Ig	Mouse Mono	WB, IP, IHC, IF, ELISA	Human, Mouse, Rat
mTOR	KDNKO 109 20657-1-AP	Rabbit Poly	WB, IP, IHC, IF, ELISA	Human, Mouse, Rat, Chicken, Cow + more
NFKB1	KDNKO 2 23576-1-AP	Rabbit Poly	WB, IHC, IF, chIP, ELISA	Human
NFKB1	66992-1-Ig	Mouse Mono	WB, IF, ELISA	Human
NFKB1.p105,p50-Specific	KDNKO 4 15506-1-AP	Rabbit Poly	WB, ELISA	Human, Rat
NFKB1	KDNKO 88 14220-1-AP	Rabbit Poly	WB, IP, IHC, IF, chIP, ELISA	Human
NF-κB p65	KDNKO 6 66535-1-Ig	Mouse Mono	WB, IHC, IF, ELISA	Human, Mouse, Rat
NF-κB p65	KDNKO 423 10745-1-AP	Rabbit Poly	WB, IP, IHC, IF, FC, chIP, ELISA	Human, Mouse, Rat, Pig, A. Sphaerocephala + more
p70(S6K)	KDNKO 60 14485-1-AP	Rabbit Poly	WB, IP, IHC, IF, FC, ELISA	Human, Mouse, Rat, Blunt Snout Bream
p70(S6K)	KDNKO 4 66638-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Mouse, Rat
PDK1	KDNKO 3 18262-1-AP	Rabbit Poly	WB, IHC, IF, ELISA	Human, Mouse, Rat
PHD1	KDNKO 4 12984-1-AP	Rabbit Poly	WB, IF, ELISA	Human, Mouse, Rat
PHD2/EGLN1	KDNKO 1 19886-1-AP	Rabbit Poly	WB, IHC, ELISA	Human, Mouse, Rat
PHD2/EGLN1	20368-1-AP	Rabbit Poly	WB, IHC, IF, ELISA	Human, Mouse
PHD2/EGLN1	66589-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Mouse, Rat, Pig
PHD3	KDNKO 3 18325-1-AP	Rabbit Poly	WB, IP, IHC, ELISA	Human
PI3K p55(gamma)	KDNKO 27035-1-AP	Rabbit Poly	WB, IHC, ELISA	Human
PI3K p110(alpha)	KDNKO 2 67071-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Mouse, Rat, Pig, Monkey
PI3K p110(beta)	KDNKO 1 21739-1-AP	Rabbit Poly	IHC, IF, ELISA	Human
PI3K p110(beta)	KDNKO 1 67121-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Rat
PI3K p85 alpha	KDNKO 28 60225-1-Ig	Mouse Mono	WB, IP, IHC, ELISA	Human, Mouse, Pig, Rat, Dog
REDD1 specific	KDNKO 162 10638-1-AP	Rabbit Poly	WB, IP, IHC, IF, chIP, ELISA	Human, Meriones Unguiculatus, Pig, Rabbit
REDD1	67059-1-Ig	Mouse Mono	WB, ELISA	Human
RPS6	KDNKO 3 14823-1-AP	Rabbit Poly	WB, IF, ELISA	Human, Mouse
RPS6	KDNKO 1 66886-1-Ig	Mouse Mono	WB, IF, ELISA	Human, Mouse, Rat
STAT3	KDNKO 10 60199-1-Ig	Mouse Mono	WB, IP, IHC, IF, ELISA	Human, Mouse, Rat
STAT3	KDNKO 13 10253-2-AP	Rabbit Poly	WB, IP, IHC, IF, FC, CoIP, ChIP, ELISA	Human, Mouse, Rat, Marmoset, Zebrafish
TLR4	KDNKO 12 66350-1-Ig	Mouse Mono	WB, IHC, ELISA	Human, Mouse, Rat, Pig
TLR4	KDNKO 148 19811-1-AP	Rabbit Poly	WB, IHC, IF, CoIP, ELISA	Human, Mouse, Rat, Calf, Hamster, Swine
TNF alpha	KDNKO 139 17590-1-AP	Rabbit Poly	WB, ELISA	Human, Mouse, Carp, Chick
TNF alpha	KDNKO 98 60291-1-Ig	Mouse Mono	WB, IHC, IF, FC, ELISA	Human, Chicken
CL488-conjugated TNF alpha	CL488-60291	Mouse Mono	IF, ELISA	Human
CL594-conjugated TNF alpha	CL594-60291	Mouse Mono	IF, ELISA	Human
VEGFA	KDNKO 7 26381-1-AP	Rabbit Poly	WB, IHC, ELISA	Human
VEGFA	66828-1-Ig	Mouse Mono	WB, IHC, IF, ELISA	Human, Mouse
VHL	24756-1-AP	Rabbit Poly	WB, IF, ELISA	Human

[KDNKO](#) the number of references to antibodies