

For Research Use Only

SYNGR2 Polyclonal antibody

Catalog Number: 10461-1-AP



Basic Information

Catalog Number:

10461-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0720

GenBank Accession Number:

BC000407

GeneID (NCBI):

9144

UNIPROT ID:

O43760

Full Name:

synaptogyrin 2

Calculated MW:

25 kDa

Observed MW:

25-27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, mouse kidney tissue, mouse lung tissue, mouse skeletal muscle tissue, rat lung tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10461-1-AP (SYNGR2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.