

For Research Use Only

TIRAP Polyclonal antibody

Catalog Number: 10497-1-AP

2 Publications



Basic Information

Catalog Number:

10497-1-AP

Concentration:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0771

GenBank Accession Number:

BC032474

GeneID (NCBI):

114609

UNIPROT ID:

P58753

Full Name:

toll-interleukin 1 receptor (TIR)
domain containing adaptor protein

Calculated MW:

256 aa, 28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Peng Li	30114378	Bosn J Basic Med Sci	WB
Xu Kedong K	23545718	Oncol Rep	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

