

For Research Use Only

Calmodulin Polyclonal antibody

Catalog Number: 10521-1-AP **2 Publications**



Basic Information

Catalog Number:

10521-1-AP

Size:

253 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0820

GenBank Accession Number:

BC006464

GeneID (NCBI):

805

UNIPROT ID:

P62158

Full Name:

calmodulin 2 (phosphorylase kinase, delta)

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Heng-Xue Wang	34484665	Toxicol Res (Camb)	WB
Jennie Sjöhamn	27643892	Protein Sci	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data