

For Research Use Only

# His-Bacteria-Protein Polyclonal antibody



Catalog Number: 10560-1-Ig **1 Publications**

## Basic Information

Catalog Number:

10560-1-Ig

Size:

773  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0050

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

30 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:2000-1:16000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : 30 kDa His-tag bacteria protein of the E.Coli(DE3),

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Ziling Fang	29547122	Elife	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data