

For Research Use Only

TH1L Polyclonal antibody

Catalog Number: 11226-1-AP **3 Publications**



Basic Information

Catalog Number:

11226-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1750

GenBank Accession Number:

BC014952

GeneID (NCBI):

51497

UNIPROT ID:

Q8IXH7

Full Name:

TH1-like (Drosophila)

Calculated MW:

66 kDa

Observed MW:

66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP, CHIP

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : A549 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Katarzyna Kaczmarek Michaels	26379089	Virology	ChIP
Matthew Gagne	31116788	PLoS Pathog	chIP
Natarajan Malini M	23884411	J Biol Chem	WB,IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 11226-1-AP (TH1L antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.