

For Research Use Only

NELL2 Polyclonal antibody

Catalog Number: 11268-1-AP **7 Publications**



Basic Information

Catalog Number: 11268-1-AP	GenBank Accession Number: BC020544	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 4753	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q99435	
Isotype: IgG	Full Name: NEL-like 2 (chicken)	
Immunogen Catalog Number: AG1793	Calculated MW: 91 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : fetal human brain tissue,
Cited Applications: WB, IF, IHC	
Species Specificity: human	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hongpan Zhang	33869274	Front Mol Biosci	IHC, WB
Fangkun Liu	30599272	Neuroscience	IF
Jianhong Fu	33416178	Oncol Rep	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 11268-1-AP (NELL2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.