

VCAM-1 Polyclonal antibody

Catalog Number: 11444-1-AP

Featured Product

83 Publications

Basic Information

Catalog Number:

11444-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1993

GenBank Accession Number:

BC017276

GeneID (NCBI):

7412

UNIPROT ID:

P19320

Full Name:

vascular cell adhesion molecule 1

Calculated MW:

739 aa, 81 kDa

Observed MW:

110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:400-1:1600

Applications

Tested Applications:

IHC, WB, ELISA

Cited Applications:

FC, IF, IHC

Species Specificity:

human

Cited Species:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: TNF alpha treated HUVEC cells,

IHC: human tonsillitis tissue, human kidney tissue, human spleen tissue, human placenta tissue, human lung cancer tissue

Background Information

Vascular cell adhesion molecule 1 (VCAM1), also known as CD106, is a 110-kDa transmembrane glycoprotein belonging to the immunoglobulin gene superfamily. VCAM1 is expressed by cytokine-activated endothelium, interacts with integrin VLA4 (α4β1) present on the surface of leukocytes, and mediates both adhesion and signal transduction. It is also expressed either constitutively or inducibly in a variety of other cell types, including vascular smooth muscle cells, differentiating skeletal muscle cells, renal and neural epithelial cells, macrophages (Kupffer cells), dendritic cells, and bone marrow stromal cells (PMID: 7507076, 11359843).

Notable Publications

Author	Pubmed ID	Journal	Application
Xinyu Zhu	36103933	Eur J Pharmacol	
Yao Ding	34492084	PLoS Pathog	IHC
Zinan Wang	34715363	Int J Biochem Cell Biol	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 4006900926



E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Untreated or TNF- α -treated HUVEC cells were subjected to SDS PAGE followed by western blot with 11444-1-AP (VCAM-1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11444-1-AP (VCAM-1 antibody) at dilution of 1:800 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).