

For Research Use Only

RAX2 Polyclonal antibody

Catalog Number: 11466-1-AP **1 Publications**



Basic Information

Catalog Number: 11466-1-AP	GenBank Accession Number: BC018709	Purification Method: Antigen affinity purification
Size: 1200 µg/ml	GeneID (NCBI): 84839	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q96153	
Isotype: IgG	Full Name: retina and anterior neural fold homeobox 2	
Immunogen Catalog Number: AG2046	Calculated MW: 184 aa, 20 kDa	
	Observed MW: 22-27 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells,
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Rui Su	32015336	Cell Death Dis	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 11466-1-AP (RAX2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.