

For Research Use Only

STAC Polyclonal antibody

Catalog Number: 11480-1-AP



Basic Information

Catalog Number:

11480-1-AP

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2059

GenBank Accession Number:

BC020221

GeneID (NCBI):

6769

UNIPROT ID:

Q99469

Full Name:

SH3 and cysteine rich domain

Calculated MW:

402 aa, 45 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, PC-3 cells, human liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 11480-1-AP (STAC antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.