

For Research Use Only

# HCK Polyclonal antibody

Catalog Number: 11600-1-AP

Featured Product

2 Publications



## Basic Information

Catalog Number:

11600-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2132

GenBank Accession Number:

BC014435

GeneID (NCBI):

3055

UNIPROT ID:

P08631

Full Name:

hemopoietic cell kinase

Calculated MW:

57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : human lung cancer tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Xiangxi Kong	33041826	Front Pharmacol	WB
Yue Zheng	35845072	Front Cardiovasc Med	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung cancer using 11600-1-AP (HCK antibody) at dilution of 1:50 (under 10x lens).