

For Research Use Only

PPP1R7 Polyclonal antibody

Catalog Number: 11664-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

11664-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2263

GenBank Accession Number:

BC000910

GeneID (NCBI):

5510

UNIPROT ID:

Q15435

Full Name:

protein phosphatase 1, regulatory
(inhibitor) subunit 7

Calculated MW:

360 aa, 40 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB: PC-3 cells, mouse brain tissue, mouse liver tissue,
HEK-293 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jingshu Tang	32642407	Acta Pharm Sin B	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 11664-1-AP (PPP1R7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.