

For Research Use Only

MOBKL1B Polyclonal antibody

Catalog Number: 11669-1-AP



Basic Information

Catalog Number:

11669-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2252

GenBank Accession Number:

BC003398

GeneID (NCBI):

55233

UNIPROT ID:

Q9H859

Full Name:

MOB1, Mps One Binder kinase
activator-like 1B (yeast)

Calculated MW:

216 aa, 25 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : K-562 cells, A375 cells, mouse liver tissue, Raji cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 11669-1-AP (MOBK1B Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.