

For Research Use Only

JWA/GTRAP3-18 Polyclonal antibody



Catalog Number: 11877-1-AP

Basic Information

Catalog Number:

11877-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2470

GenBank Accession Number:

BC020797

GeneID (NCBI):

10550

UNIPROT ID:

O75915

Full Name:

ADP-ribosylation-like factor 6
interacting protein 5

Calculated MW:

188 aa, 22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

This protein, ARL6IP5 or GTRAP3-18, localizes to the cell membrane and cytoplasm, and specifically interacts with the carboxy-terminal intracellular domain of EAAC1, functioning in EAAC1-mediated glutamate transport. The monomer tended to form multimers including the dimer(50-60kd).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data