

For Research Use Only

TRA2A Polyclonal antibody

Catalog Number: 12079-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

12079-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2760

GenBank Accession Number:

BC017094

GeneID (NCBI):

29896

UNIPROT ID:

Q13595

Full Name:

transformer 2 alpha homolog
(Drosophila)

Calculated MW:

282 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : A549 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Wenpeng Cao	35635094	Bioengineered	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 12079-1-AP (TRA2A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.