

For Research Use Only

# BFL-1/A1 Polyclonal antibody

Catalog Number: 12223-1-AP **2 Publications**



## Basic Information

<b>Catalog Number:</b> 12223-1-AP	<b>GenBank Accession Number:</b> BC016281	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 270 µg/ml	<b>GeneID (NCBI):</b> 597	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q16548	
<b>Isotype:</b> IgG	<b>Full Name:</b> BCL2-related protein A1	
<b>Immunogen Catalog Number:</b> AG2844	<b>Calculated MW:</b> 175 aa, 20 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB

**Species Specificity:**  
human

**Cited Species:**  
human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yan-Wei Hu	26201458	Apoptosis	WB
Ming-Da Wang	26698170	Hepatology	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data