

For Research Use Only

# FHL5 Polyclonal antibody

Catalog Number: 12650-2-AP



## Basic Information

Catalog Number:

12650-2-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3364

GenBank Accession Number:

BC021723

GeneID (NCBI):

9457

UNIPROT ID:

Q5TD97

Full Name:

four and a half LIM domains 5

Calculated MW:

284 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse heart tissue, human heart tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 12650-2-AP (FHL5 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.