

For Research Use Only

# SERPIND1 Polyclonal antibody

Catalog Number: 12741-1-AP **2 Publications**



## Basic Information

Catalog Number:

12741-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3637

GenBank Accession Number:

BC035028

GeneID (NCBI):

3053

UNIPROT ID:

P05546

Full Name:

serpin peptidase inhibitor, clade D  
(heparin cofactor), member 1

Calculated MW:

499 aa, 57 kDa

Observed MW:

57-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : human plasma tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jo-Chuan Liu	33953300	Sci Rep	WB
Liancheng Zhu	27430660	Oncol Rep	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



human plasma tissue were subjected to SDS PAGE followed by western blot with 12741-1-AP (SERPIND1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.