

For Research Use Only

ELP4 Polyclonal antibody

Catalog Number: 12746-1-AP



Basic Information

Catalog Number:

12746-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3707

GenBank Accession Number:

BC012514

GeneID (NCBI):

26610

UNIPROT ID:

Q96EB1

Full Name:

elongation protein 4 homolog (S.

cerevisiae)

Calculated MW:

424 aa, 46 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IHC 1:20-1:200

Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: human cerebellum tissue, human brain tissue

IHC: human endometrial cancer tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



human cerebellum tissue were subjected to SDS PAGE followed by western blot with 12746-1-AP (ELP4 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human endometrial cancer using 12746-1-AP (ELP4 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human endometrial cancer using 12746-1-AP (ELP4 antibody) at dilution of 1:50 (under 40x lens).