

For Research Use Only

SRGAP3 Polyclonal antibody

Catalog Number: 12823-1-AP **1 Publications**



Basic Information

Catalog Number: 12823-1-AP	GenBank Accession Number: BC039300	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 9901	
Source: Rabbit	UNIPROT ID: O43295	
Isotype: IgG	Full Name: SLIT-ROBO Rho GTPase activating protein 3	
Immunogen Catalog Number: AG3516	Calculated MW: 1099 aa, 125 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IHC

Species Specificity:
human, mouse, rat

Cited Species:
rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Anjie Lu	25206756	Neural Regen Res	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data