

For Research Use Only

EME1 Polyclonal antibody

Catalog Number: 12975-1-AP



Basic Information

Catalog Number:

12975-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3580

GenBank Accession Number:

BC016470

GeneID (NCBI):

146956

UNIPROT ID:

Q96AY2

Full Name:

essential meiotic endonuclease 1

homolog 1 (S. pombe)

Calculated MW:

583 aa, 65 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: A375 cells, human liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A375 cells were subjected to SDS PAGE followed by western blot with 12975-1-AP (EME1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.