

For Research Use Only

PASD1 Polyclonal antibody

Catalog Number:12976-1-AP



Basic Information

Catalog Number: 12976-1-AP	GenBank Accession Number: BC040301	Purification Method: Antigen affinity purification
Size: 750 ug/ml	GeneID (NCBI): 139135	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	UNIPROT ID: Q8IV76	
Isotype: IgG	Full Name: PAS domain containing 1	
Immunogen Catalog Number: AG3587	Calculated MW: 773 aa, 87 kDa	
	Observed MW: 72-87 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : NCI-H1299 cells, U2OS cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12976-1-AP (PASD1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.