

For Research Use Only

DBF4 Polyclonal antibody

Catalog Number: 13076-1-AP



Basic Information

Catalog Number:

13076-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3771

GenBank Accession Number:

BC036045

GeneID (NCBI):

10926

UNIPROT ID:

Q9UBU7

Full Name:

DBF4 homolog (S. cerevisiae)

Calculated MW:

77 kDa

Observed MW:

70-80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : NIH/3T3 cells, Jurkat cells, PC-3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 13076-1-AP (DBF4 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.