

For Research Use Only

NUP93 Polyclonal antibody

Catalog Number: 13077-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

13077-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3735

GenBank Accession Number:

BC034346

GeneID (NCBI):

9688

UNIPROT ID:

Q8N1F7

Full Name:

nucleoporin 93kDa

Calculated MW:

819 aa, 93 kDa

Observed MW:

85-93 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, rat, mouse

Positive Controls:

WB : mouse testis tissue, A431 cells, HEK-293 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhixiong Liu	30876848	Neuron	WB
Yuri Furusawa	30087672	Front Microbiol	WB
Bin Chen	35152263	Oncogene	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13077-1-AP (NUP93 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.