

For Research Use Only

RBMS2 Polyclonal antibody

Catalog Number: 13166-1-AP **1 Publications**



Basic Information

Catalog Number:

13166-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3818

GenBank Accession Number:

BC027863

GeneID (NCBI):

5939

UNIPROT ID:

Q15434

Full Name:

RNA binding motif, single stranded interacting protein 2

Calculated MW:

407 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : mouse skin tissue, HeLa cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
N Pospiech	29434205	Sci Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 13166-1-AP (RBMS2 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.