

For Research Use Only

# CD32 Polyclonal antibody

Catalog Number: 13173-1-AP



## Basic Information

Catalog Number:

13173-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3828

GenBank Accession Number:

BC019931

GeneID (NCBI):

2212

UNIPROT ID:

P12318

Full Name:

Fc fragment of IgG, low affinity IIa, receptor (CD32)

Calculated MW:

316 aa, 35 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data