

For Research Use Only

ALOX15 Polyclonal antibody

Catalog Number: 13286-1-AP



Basic Information

Catalog Number:

13286-1-AP

Size:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4040

GenBank Accession Number:

BC029032

GeneID (NCBI):

246

UNIPROT ID:

P16050

Full Name:

arachidonate 15-lipoxygenase

Calculated MW:

662 aa, 75 kDa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse colon tissue, rat colon tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13286-1-AP (ALOX15 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.