

For Research Use Only

TRAM2 Polyclonal antibody

Catalog Number: 13311-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

13311-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4144

GenBank Accession Number:

BC028121

GeneID (NCBI):

9697

UNIPROT ID:

Q15035

Full Name:

translocation associated membrane protein 2

Calculated MW:

370 aa, 43 kDa

Observed MW:

43 kDa, 66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HepG2 cells, Sp2/O cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Fangxing Xu	36618348	Front Immunol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 13311-1-AP (TRAM2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.