

For Research Use Only

IKZF3 Polyclonal antibody

Catalog Number: 13561-1-AP **2 Publications**



Basic Information

Catalog Number:

13561-1-AP

Size:

367 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4434

GenBank Accession Number:

BC032707

GeneID (NCBI):

22806

UNIPROT ID:

Q9UKT9

Full Name:

IKAROS family zinc finger 3 (Aiolos)

Calculated MW:

509 aa, 58 kDa

Observed MW:

58-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB: Jurkat cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Konstantinos Dimopoulos	29718735	Leuk Lymphoma	IHC
Liang Zhou	33392195	Front Cell Dev Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 13561-1-AP (IKZF3 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.