

For Research Use Only

Aiolos/IKZF3 Polyclonal antibody

Catalog Number: 13561-1-AP **3 Publications**



Basic Information

Catalog Number: 13561-1-AP	GenBank Accession Number: BC032707	Purification Method: Antigen affinity purification
Size: 367 µg/ml	GeneID (NCBI): 22806	Recommended Dilutions: WB 1:500-1:2400
Source: Rabbit	UNIPROT ID: Q9UKT9	
Isotype: IgG	Full Name: IKAROS family zinc finger 3 (Aiolos)	
Immunogen Catalog Number: AG4434	Calculated MW: 509 aa, 58 kDa	
	Observed MW: 58-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells,
Cited Applications: WB, IHC	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Konstantinos Dimopoulos	29718735	Leuk Lymphoma	IHC
Liang Zhou	33392195	Front Cell Dev Biol	WB
Busong Wang	39136771	Cell Mol Life Sci	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 13561-1-AP (IKZF3 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.