

For Research Use Only

# ICK Polyclonal antibody

Catalog Number: 13611-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 13611-1-AP	<b>GenBank Accession Number:</b> BC035807	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 130 µg/ml	<b>GeneID (NCBI):</b> 22858	<b>Recommended Dilutions:</b> WB 1:500-1:2000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UPZ9	
<b>Isotype:</b> IgG	<b>Full Name:</b> intestinal cell (MAK-like) kinase	
<b>Immunogen Catalog Number:</b> AG4511	<b>Calculated MW:</b> 632 aa, 71 kDa	
	<b>Observed MW:</b> 71 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : PC-3 cells,
<b>Cited Applications:</b> IF	
<b>Species Specificity:</b> human, mouse, rat	
<b>Cited Species:</b> mouse	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
I-Ling Lu	29950674	Nat Commun	IF

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 13611-1-AP (ICK antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.